

ConvaQuip Ind., Inc.

Manual for Model 724A Bedside Commode



Capacity 750 lbs.

Model 724A:

CAUTION: Make sure toilet seat ring is firmly in the down position and resting against the top of the frame in order to prevent pinching skin. Once seat ring is in place, slowly lower yourself into a sitting position. Keep hands and fingers on the armrests while lowering into position.

CAUTION: Position the bedside commode against a wall prior to use. The tubular back needs to be touching the wall when entering or exiting the bedside commode to prevent tipping.

Cleaning Instructions:

- * For infection control, clean all parts with one part bleach and 9 parts water. Rinse and allow to air dry or wipe with a clean, dry cloth.
- * Clean rubber tips with alcohol to remove oil, grime and dirt
- * Pail may be steam cleaned or use Beaucoup disinfectant (made by EcoLab) or equivalent.

ConvaQuip Ind., Inc.

Assembly

Model 724A
Bedside Commode



Set top half on bottom half.
Slide down over the tubing until the
tubing will rest on all four snap buttons.



Close up of tubing resting on snap button.

Press in the snap button on each leg and barley push down the tubing, just enough for the button to disappear. Now slowly and in small increments move from leg to leg pushing the tubing down.

If you push too far on any leg, it will bind up and you will have to even them out again by pulling it up and starting over. Once all four tubes have bottomed out, the button will reappear and snap into place. When all buttons are protruding through the button hole, the commode is ready to use.

